

Quality Notification

QN-RIS-DiagnosticsInsights-2025-004

RIS/Roche Information Solutions
Post-Market Quality
Version 1
11-Jul-2025

"eGFR" test received from the NovaNet StatStrip instrument is mapped to the "CrCl" test and the "CrCl" result test received from the instrument is not processed in the navify POC Operations application

Product	GMMI
cobas infinity POC core license - 2.x	09159053001
navify POC Operations - 3.x	10170280001
Instrument / System	navify® POC Operations
Components	Software Modification Tracking
Target Group	Service Application
Details	Feedback required by Affiliate Field Implementation Mandatory Next site visit / Implementation Date Affiliate action
Implementation Date	31-Dec-2025
Feedback Due Date	30-Oct-2025
Escalated Case Number	CN-1045635
CAPA/NC/SCAR	N/A

"eGFR" test received from the NovaNet StatStrip instrument is mapped to the "CrCl" test and the "CrCl" result test received from the instrument is not processed in the navify POC Operations application

Modification Label	S_ZH2025/004Q Biomedical NovaNet driver
Duration Estimate	40 Minutes

Executive Summary

This Quality Notification informs about an issue encountered regarding the “eGFR” test received from the NovaNet StatStrip instrument, which is mapped to the “CrCl” test in the navify POC Operations application and the “CrCl” test received from the instrument is not processed.

This issue has been found in the Nova Biomedical NovaNet version 1.1.1 driver used for connecting NovaNet StatStrip instrument to the cobas infinity POC/navify POC Operations application and has been fixed in the already released version Nova Biomedical NovaNet version 1.2.0.

Units Concerned

Nova Biomedical NovaNet version 1.1.1 driver.

Subject

An allegation was received indicating that, despite the correct configuration of the “eGFR” test, its results never appeared in the navify POC Operations application. Meanwhile, other tests, such as Creatinine and Clearance “CrCl”, were displayed correctly in the application.

During the investigation of the complaint case CN-1045635, it was found that based on the StatStrip instrument Host Interface Manual from NovaNet, the test code “CREGFR” sent by the device should be assigned to the “eGFR” test name in the application. However, it was found that for the test code “CREGFR” the transformation was done to the “CrCl” test name instead of the “eGFR” test name, as is described in the Host Interface Manual of the instrument. Furthermore, it was noticed that when receiving the “CRESS” test code from the device, it was correctly mapped by the driver to the “Creat” test name in the application.

For this allegation about the “eGFR” not being shown in the navify POC Operations application and the Creatinine and Clearance (CrCl) working fine instead, it could be found that this was caused by two issues in the driver:

"eGFR" test received from the NovaNet StatStrip instrument is mapped to the "CrCl" test and the "CrCl" result test received from the instrument is not processed in the navify POC Operations application

- When the test code "CREGFR" is received from the NovaNet instrument, instead of being mapped to the "eGFR" Test name, the driver is mapping it to the "CrCl" test.
- The "CRECRCL" test code received from the NovaNet instrument is not mapped to any test in the application.

Hence, it could be seen that the "eGFR" test was present in the application but using the "CrCl" Test name instead and that the "CRECRCL" which should be mapped to the "CrCl" was not visible in the application.

Root Cause

The driver is incorrectly mapping the "CREGFR" received from the NovaNet instrument (mapping it to the "CrCl" instead of to the "eGFR" Test name). In parallel, the "CRECRCL" test code received from the NovaNet instrument is not mapped to any test in the application (when it should be mapped to the "CrCl").

Risk Assessment

The Risk has been considered as "Moderate" due to the incorrect mapping of a test that leads to instrument results being applied to a different test.

Frequency of Occurrence

The Nova Biomedical NovaNet 1.1.1 was released in Q1 2024; since then, no related cases have been reported until this first allegation (CN-1045635).

Detectability

It was highly detectable that the test was not present in the application, but not easily detectable that it was mapped to a different test, as this other test (CrCl) received from the instrument was itself not mapped. However, it could also be observed that the Units differ. For the "eGFR" test, the units used are "mL/min/1.73 m²" and in the case of the Creatinine Clearance "CrCl", the units used are "mL/min", which makes the detectability higher.

Severity

No diagnostic test results are affected and a medical risk to patients and users can be excluded. Therefore, a Health Hazard Evaluation is not required.

"eGFR" test received from the NovaNet StatStrip instrument is mapped to the "CrCl" test and the "CrCl" result test received from the instrument is not processed in the navify POC Operations application

Important Information

The Nova Biomedical NovaNet version 1.1.1 driver needs to be updated to the Nova Biomedical NovaNet version 1.2.0 driver.

Instructions

To fix the issue, upgrade the Nova Biomedical NovaNet version 1.1.1 driver to the Nova Biomedical NovaNet version 1.2.0 driver.

Actions

Internal Actions	The issues have been tracked with the following defects: POCC-839 and POCC-841
Actions for Affiliates	- To upgrade the Nova Biomedical NovaNet version 1.1.1 driver to the Nova Biomedical NovaNet version 1.2.0 driver - Please submit the QN feedback via the Smarteeva IT tool by the due date - Please use the provided modification label (change order) to document the installation of the latest driver. Note that no work orders will be automatically generated as this driver is used at a limited number of customer sites

Contact

Please do not hesitate to contact us in case of questions regarding the information provided.

Best Regards,

Author

Marta Güell Pasero
Senior Engineer PMQ

Reviewer

Yerko Yanzich
Principal Engineer PMQ