

Support Notification

SN-RIS-DiagnosticsInsights-2025-096

RIS/Roche Information Solutions

Post-Market Quality

Version 1

17-Jul-2025

navify POC Operations and cobas infinity POC statement about InterSystems Iris Security Notice SVR2025A

Product	GMMI
navify POC Operations - 3.x	10170280001
cobas infinity POC core license	09159053001
Instrument / System	navify® POC Operations cobas® infinity POC Add-on
Components	Software
Target Group	Service
Details	Internal information Risk Statement included
Escalated Case Number	CN-1059105
CAPA/NC/SCAR	NPOCO-23697

Executive Summary

The recent announcement of vulnerabilities (SVR2025A) by InterSystems outlines a natively **high criticality level** due to potential security implications on their platform. However, after a thorough assessment of these vulnerabilities concerning our products, it has been determined that their effective criticality for **navify** POC Operations and **cobas infinity** POC is **low**. This is primarily due to how these solutions are deployed on secured hospital infrastructure, isolated from public networks.

navify POC Operations and cobas infinity POC statement about InterSystems Iris Security Notice SVR2025A

InterSystems has shared the announcement with its application partners and direct end-user customers. A broader public announcement has been postponed to provide InterSystems customers sufficient time to develop and implement remediation plans where necessary.

This notification is intended for affiliates that require local authorities to issue mitigations to vendors or provide official statements when the disclosed vulnerabilities are considered critical in severity and have become public.

Units Concerned

Servers running **cobas infinity** POC.
Servers running **navify** POC Operations.

Subject

The recent announcement of vulnerabilities (SVR2025A) by InterSystems outlines a natively **high criticality level** due to potential security implications on their platform.

SVR2025A-1:

Class methods in web applications may execute actions without proper privilege checking.
CVSS 3.1 Score/Severity: **9.4 Critical**.

SVR2025A-2:

Malicious links to web applications, if clicked, may execute arbitrary code.
CVSS 3.1 Score/Severity: **8.3 High**.

However, upon careful evaluation of these vulnerabilities in the context of our products conducted by our Product Security and Privacy Office (PSPO) team in collaboration with the R&D team, it has been determined that their effective criticality for **navify** POC Operations and **cobas infinity** POC is **Low** for the following reasons:

Mitigation by Deployment and Use Context: The sole external attack vector is limited to the use of a public IIS web server, but this does not directly expose the InterSystems IRIS web applications connected to our products. Furthermore, these solutions are deployed within secure hospital infrastructure and are not exposed to public networks. Exploitation of these vulnerabilities would require an attacker to gain initial operating system level access to the server itself, which represents an entirely different class of security breach. In such an event, the integrity of the entire server (including the database) would already be compromised, rendering the SVR2025A vulnerabilities a secondary concern.

navify POC Operations and cobas infinity POC statement about InterSystems Iris Security Notice SVR2025A

Protection Measures in Place: Our products already include key user assistance recommendations to enhance overall server and system security, explained in the User Assistance:

- Make sure other computers and services on the network (for example, the LIS, archiving share, backup share, or service) are properly secured and protected against malicious software and unauthorized access.
- Restrict physical access to the system and server and all attached IT infrastructure (computer, cables, network equipment, and so on).
- Ensure that all device components maintaining private data are physically secure (i.e. require a physical key to be accessed or removed).

These measures align with industry-standard best practices to minimize access vectors to critical systems.

Planned fix Implementation: The project has decided to natively address these vulnerabilities by incorporating patches provided by InterSystems into the upcoming **navify** POC Operations version 3.3.0, which is planned for release in the fourth quarter (Q4) of this year.

Risk Statement

Based on the current knowledge of the disclosed vulnerability and the thorough evaluation conducted by our PSPO team in collaboration with the R&D team, a risk to patients and users can be excluded.

Parts Required

N/A

Further Information

N/A

Instructions

navify POC Operations and cobas infinity POC statement about InterSystems Iris Security Notice SVR2025A

It is suggested to use the information within this document in case it is required to present an official statement about vulnerabilities classified as critical to your local authorities.

Actions

Internal Actions

The vulnerabilities reported by InterSystems in the SVR2025A Alert are being addressed in the defect NPOCO-23697, currently assigned to be fixed in **navify** POC Operation version 3.3.0.

Actions for Affiliates

It is recommended to use this document as a reference if you need to present an official statement about vulnerabilities classified as critical to your local authorities.

Contact

Please do not hesitate to contact us in case of questions regarding the information provided.

Best Regards,

Author

Yerko Yanzich

Principal Engineer PMQ

Reviewer

Jörg-Dieter Schlederer

GCS RIS Lab Insights Network TPM

Reviewer

Jakub Winiarz

Business Product Manager - Point of Care software