

Support Notification

SN-RIS-DiagnosticsInsights-2025-094

RIS/Roche Information Solutions

Global Customer Support

Version 2

08-Aug-2025

navify POC Operations 3.2 released globally

Product	GMMI
cobas infinity POC core license - 2.x	09159053001
navify POC Operations - 3.x	10170280001
Instrument / System	navify ® POC Operations cobas ® infinity POC Add-on
Components	Modification Tracking Documentation Software
Target Group	Service
Details	Internal information Affiliate action Risk Statement included
SW Version Tracking	POC v3.2
CAPA/NC/SCAR	N/A
Attachments	Software Release Notes
Modification Label	S_ZH2025/094S navify POC Operations 3.2 S_H22025/094S navify POC Operations 3.2
Duration Estimate	30 Minutes

navify POC Operations 3.2 released globally

Change History

Version	Date	Description
1	14-Jul-2025	Original Version
2	08-Aug-2025	New language pack installer

Executive Summary

navify POC Operations V3.2 refines and enhances the functionality provided by **navify** POC Operations V3.1. It also serves as the functional replacement for the **cobas IT** 1000 frontend, frees the **cobas POC** IT bundle of products from the dependency on Microsoft Silverlight as it is built on HTML 5. It incorporates frontend functionality as well as the backend previously provided by the **cobas IT** 1000 application backend.

Units Concerned

navify POC Operations V3.2 can be installed where there is currently no version of **navify** POC Operations, **cobas infinity** POC or **cobas IT** 1000 installed or as a direct upgrade to released versions higher or equal to **cobas IT** 1000 V2.10.00 and **cobas infinity** POC 1.0.1, older versions need to be updated to **cobas IT** 1000 V2.10.00 or newer first. Thus far, **navify** POC Operations V3.2 is available in English (out of box), language packs for further languages will be released and published as soon as possible.

GCS Product Support Desk (PSD)

<https://docs.google.com/spreadsheets/d/1Polck_UBsJstg1cH3dfaCHdxAFL-OQznPge5vrVS0J8/edit#gid=0>

Channel	Details
Hotline Phone Number (Emergency)	+41 41 799 41 41, 142 (POC IT)
REXIS (Support requests, Complaints and Inquiries)	Click 'New Global Support Ticket' in Escalated Issues section
Mailbox (License and Training)	global.gcs_rpd-pocit@roche.com
Google Chat (Questions)	global.gcs_rpd-pocit@roche.com

Note: PSD coverage is available Monday – Friday, 09:00 - 17:00 (CET) Central European Time Zone.

navify POC Operations 3.2 released globally

Complaint Handling

Log complaints and inquiries on the following parts as applicable:

Global Material Number	Description
09159053001	cobas infinity POC core license - 2.x
10170280001	navify POC Operations - 3.x

Subject

For a list of

- new features
- changes and solved issues
- released items
- limitations
- known open issues and workarounds
- constraints see the attached Software Release Notes documents or a summary of them at <https://rdkm.roche.com/backgrounds/484166>.

Release notes are for Roche internal use only and must not be shared with the end-customer directly.

navify POC Operations 3.2 is supported on the Google Chrome and Microsoft Edge browsers for the frontend on PCs and on iPadOS and Android for mobile platforms. It is supported on Windows Server 2016, Windows Server 2019 and Windows Server 2022 for the backend, Windows Server 2012 R2 is no longer a supported platform. Existing Caché license keys for **cobas infinity** POC and **cobas IT 1000** must be upgraded to InterSystems IRIS for Health 2024.1 compatible licenses in order to install **navify** POC Operations V3.2.

The adapter used by the **cobas infinity** POC mobile and **cobas infinity** POC tablet mobile apps is not compatible with **navify** POC Operations and must not be installed on **navify** POC Operations servers. Consequently, the use of the **cobas infinity** POC mobile and **cobas infinity** POC tablet mobile apps is not supported with **navify** POC Operations.

cobas b 123 and **cobas b 221** instruments connected to **navify** POC Operations V3.2 in conjunction with **cobas infinity** POC remote access cannot be connected to the **cobas bge** link software in parallel. The **cobas bge** link server component cannot be active on the server hosting the **cobas infinity** POC remote access instrument connectivity if **cobas b 221** instruments are connected to **cobas infinity** POC remote access.

navify POC Operations 3.2 released globally

Risk Statement

No patient or diagnostic test results are affected and a medical risk to patients and users can be excluded.

Training

As **navify** POC Operations 3.2 technologically is a continuation of **cobas infinity** POC V2, **cobas infinity** POC V2 content remains valid for **navify** POC Operations 3.2, the following Coached Distance Learning training courses are already available in Cornerstone or will be newly available in Cornerstone shortly:

- **cobas infinity** POC V2.0 Operator Training (standardized)
- **cobas infinity** POC V2.0.1 Delta Operator Training (standardized)
(this training course also covers the differences between V2.0.0 and V2.0.1 in terms of service topics which in total are not extensive enough to warrant a separate training course)
- **cobas infinity** POC V2.1 Delta Operator Training (standardized)
(this training course also covers the differences between V2.0.1 and V2.1 in terms of service topics which in total are not extensive enough to warrant a separate training course)
- **cobas infinity** POC V2.2.x Delta Operator Training (standardized)
(this training course also covers the differences between V2.1 and V2.2.2 in terms of service topics which in total are not extensive enough to warrant a separate training course)
- **cobas infinity** POC V2.3 Operator Training Update (standardized)
(this training course also covers the differences between V2.2.2 and V2.3 in terms of service topics which in total are not extensive enough to warrant a separate training course)
- **navify** POC Operations V3.0 Operator Training Update (standardized)
(this training course also covers the differences between V2.3 and V3.0 in terms of service topics which in total are not extensive enough to warrant a separate training course)
- **navify** POC Operations V3.1 Operator Training Update (standardized)
(this training course will also cover the differences between V3.0 and V3.1 in terms of service topics which in total are not extensive enough to warrant a separate training course)
- **navify** POC Operations V3.2 Operator Training Update (standardized)
(this training course will also cover the differences between V3.0 and V3.1 in terms of service topics which in total are not extensive enough to warrant a separate training course)
- **cobas infinity** POC V2.0 Service Training (standardized)

navify POC Operations 3.2 released globally

Parts Required

Software

- **navify** POC Operations V3.2 (backend/foundation)
Installer Roche.navify.POC.Operations.Foundation_3_2_0_849.exe is published in
<https://drive.google.com/drive/folders/1WYf1H8WPJ_DhgBB0dTBPExyhgZZ-gF_J>
- **navify** POC Operations V3.2 (frontend)
Installer Roche.navify.POC.Operations.Installer.3.2.0.2530.exe is published in
<https://drive.google.com/drive/folders/1WYf1H8WPJ_DhgBB0dTBPExyhgZZ-gF_J>

Additional Software

- **Language Pack Installer**
Installer **solutionsuite-lpinstaller_2.2.0.zip** is published in
<https://drive.google.com/drive/folders/1k3TCytLTf-xyq9dJ3rY5yyn_xia_jWUK>
- **Language Packs**
The Language Packs are/will be published in
<<https://drive.google.com/drive/folders/1puzFcUjjkgzldjszkbDY4UhizkYVz6R>>

Documentation

The following documents are available via PIM.web:

- **navify** POC Operations 3.2 User Guide/User Assistance 6.2 EN (pdf)
 - **navify** POC Operations Host Interface Manual– software version 3.0– publication version 1.4 (pdf)

The following document will be available shortly via PIM.web:

- **navify** POC Operations V3.0 Host Interface Configuration Manual (pdf)

Until its publication, please refer to the respective document for the **cobas IT 1000** application.



navify POC Operations 3.2 released globally

Instructions

See the RDK Service Documentation (link [here](#)) for installation and update instructions.

Actions

Internal Actions	N/A
Actions for Affiliates	Please use the modification labels provided in this SN to report the installation through a Service Activity.

Contact

Please do not hesitate to contact us in case of questions regarding the information provided.

Best Regards,

Author

Javier Arocha Herrera

GCS RIS Support & Training Specialist

Reviewer

Jörg-Dieter Schleder

GCS RIS Lab Insights Network TPM

Affected Serial Numbers

All Units Concerned